

LabJack UE9 User's Guide Revision History

V1.00 released 4/1/2005

V1.01 released 4/11/2005

Section 2.10.2 – Minor changes.

Section 2.12.1 – New section describing the CB37 terminal board.

Section 3.2.1 – New section describing external triggering for stream mode.

Sections 4.3.13 - 4.3.15 – New sections declaring low level watchdog and reset functions.

V1.01b released 8/17/2005

Section 2.9 – Added information to the end of section.

Section 4.3.15 – Fixed soft/hard reset bits documentation.